

Observing Spin Waves

PHAROS

Benjamin G. Ueland
David Schmitter
Zhixian Zhou
Hyunjeong Kim
Ju-Hyun Park
Wei Tian
Abdel Alsmadi
Stephen Wilson
Victor Fanelli

Acknowledgements:

Frans Trouw

Markus Hehlen


Ray Osborn

The LANSCE Neutron
Scattering Winter School

Outline


- What we measured
- Why we measured it
- How we measured it
- Our results

Crystal structure of HoFe_2



HoFe_2 forms
Cubic Laves
Phase compound
with the easy axis
of magnetization
along the $[100]$
direction

Ferrimagnetic State



Magnetization of HoFe₂




FIG. 2. Hysteresis loop of TbFe₂ at 4.2°K.




Fig. 3 Spontaneous magnetic moment of single crystal HoFe₂ as a function of temperature.

$T_c = 606$ K for Single crystal

Magnon dispersion relations

$$\hbar \omega_q^{ferro} = S [J(0) - J(q)] \approx D q^2$$

In cubic symmetry
 $D = 2JSa^2$

$$\hbar \omega_q^{antiferro} = S \sqrt{J(0)^2 - J(q)^2}$$


where

$$J(q) = \sum_{i,j} J_{ij} \exp \left\{ i \vec{q} \cdot (\vec{R}_i - \vec{R}_j) \right\}$$



- Tb goes into the Ho sites
- Curie Temperature, $T_c = 610 \text{ K}$
- Inelastic neutron scattering \rightarrow Magnon Dispersion Relation \rightarrow direct information on the Exchange Interactions

Pharos Spectrometer



$\Delta E/E = 2\% \text{ to } 4\%$

Sample Size: 5 cm to 7.5 cm

Number of Detectors: 376 He³ with
40 pixels each


2-Theta : -10° to 145°

Energy Selection and Resolution

- Fixed Incident Energy
- Scan over Energy Transfers
- Choosing Incident Energy

$$TOF(\mu s / m) = \sqrt{\frac{5.2276 \cdot 10^6}{E(meV)}}$$

- Phasing of Fermi and T_0 choppers
 - For 25meV Time to Fermi Chopper = 8231 μ s
 - Tune the Frequency of Fermi Chopper to select the Energy Resolution




Sample Holder and Alignment

- Experimental Conditions:
 - Evacuated sample environment (reduces air scatter)
 - Versatile sample environment
 - Temperature, pressure, and magnetic field
- Sample Mounting and Alignment:
 - Aligned in (1,1,1) with white beam (Fermi Chopper removed)
 - Aluminum sample holder
 - Room temperature experiment
 - 1 day scan




One Slice in Time




- 40 vertical pixels in each detector bank
- Debye-Scherrer Cones from the aluminum holder intersect at straight lines, at certain angles
- Single crystal sample gives spots

Diffraction Pattern




Inelastic vs. Elastic



- Clearly see energy transfer
- Integrated over the 40 pixels in each detector bank
- Imposing a little deviation on 2σ
- Not taking into account that we are interested in a spot

Ahh...Dispersion



Ahh...Dispersion





FIG. 1. Magnon dispersion relation for $\text{Ho}_{0.38}\text{Tb}_{0.62}\text{Fe}_2$ measured at room temperature.


—R. M. Nicklow et al. *PRL* **36**, 532 (1976).

•Excellent Agreement!!!

Conclusions

- Inelastic Neutron Scattering is an excellent probe of fundamental excitations in matter
 - PHAROS offers excellent energy resolution over a wide energy range.
- Our results agree quite well with those previously reported
 - R. M. Nicklow et al. *PRL* **36**, 532 (1976).
 - J.J. Rhyne et al. *Physica* **86-88B**, 149 (1977).

THANK YOU



- The minimum Q to see Energy Gain is at Q_{elastic} .